

In Igor Pro 9 and later, you can save HDF5 packed experiment files in compressed form. This section provides a brief introduction to HDF5 compression as it relates to saving experiment files. For details on compression, see **HDF5 Compression** on page II-213.

Compression can reduce file size but also takes more time. Whether compression is worthwhile for you depends on the nature of your data and how you trade off file size versus time. If your experiments contain lots of relatively small waves with lots of noise, compression will save little disk space. If your experiments contain large waves with large sections containing one value (for example, an image that is mostly black), compression can save considerable disk space. The only way to tell if compression is worthwhile for you is to experiment with it.

To save an HDF5 packed experiment with compression via the File menu you need to enable HDF5 default compression and then provide three parameters via the Experiment section of the Miscellaneous Settings dialog. The parameters are:

- The minimum size a wave must be before compression is used
- A compression level from 0 (no compression) to 9 (max)
- Whether you want to enable shuffle (an optional process before compression)

Here are instructions for testing HDF5 compression with your experiment files. These instructions use the SaveExperiment operation via the TestHDF5ExperimentCompression function, not File->Save Experiment, and consequently are independent of the default compression settings in the Miscellaneous dialog.

1. Activate the "Test HDF5 Experiment Compression.ipf" procedure file as a global procedure file and restart Igor. The file is in WaveMetrics Procedures/File Input Output in your Igor Pro folder.
2. If you don't know how to do this, see **Activating WaveMetrics Procedure Files** on page II-33.
3. Open a typical experiment file.
4. Execute this:

```
TestHDF5ExperimentCompression(10000, 2, 0)
```

10000 is the minimum number of elements in a wave to be compressed. Smaller waves are saved uncompressed.

2 is the zip compression level.

0 means shuffle is off which is usually what you want.

The TestHDF5ExperimentCompression command saves copies of the current experiment in uncompressed form and in compressed form using the specified parameters and prints a message in the command window history area showing the effect of compression on the time to save the experiment and on the file size. The command deletes the copies so there is no junk left over.

Try different parameters for the minimum wave size and zip compression level to see how they affect save time and compression ratio. In most cases, higher zip compression levels provide small increases in compression and are not worth the time required.

Unless you find significant compression ratios, we recommend that you eschew compression. Disk space is abundant and time is precious.

By default, compression is disabled when you use the File menu to save an HDF5 packed experiment file. If you decide that you want to use compression, see **HDF5 Default Compression** on page II-214 to learn how to enable compression via the Miscellaneous Settings dialog.

### Saving as an Unpacked Experiment File

In the unpacked format, an experiment is saved as an **experiment file** and an **experiment folder**. The file contains instructions that Igor uses to recreate the experiment while the folder contains files from which Igor loads data. The experiment folder is also called the **home folder**.

## Chapter II-3 — Experiments, Files and Folders

The main utility of this format is that it is faster for experiments that contain very large numbers of waves (thousands or more). However the unpacked format is more fragile and thus is not recommended for routine use.


To save a new experiment in the unpacked format, choose Save Experiment from the File menu. At the bottom of the resulting Save File dialog, choose Unpacked Experiment Files from the popup menu. When you click Save, Igor writes the unpacked experiment file which as a ".uxp" extension.

Igor then automatically generates the experiment folder name by appending "Folder" or the Japanese equivalent, to the experiment file name. It then creates the unpacked experiment folder without further interaction. For example, if you enter "Test.uxp" as the unpacked experiment file name, Igor automatically uses "Test Folder", or the Japanese equivalent, as the unpacked experiment folder name.

If a folder named "Test Folder" already exists then Igor displays an alert asking if you want to reuse the folder for the unpacked experiment.

If the automatic generation of the unpacked experiment folder name causes a problem for you then you can save an experiment with the names of your choice using the **SaveExperiment** /F operation.

This illustration shows the icons used with an unpacked experiment and explains where things are stored.



You normally have no need to deal with the files inside the experiment folder. Igor automatically writes them when you save an experiment and reads them when you open an experiment.

If the experiment includes data folders (see Chapter II-8, **Data Folders**) other than the root data folder, then Igor will create one subfolder in the experiment folder for each data folder in the experiment. The experiment shown in the illustration above contains no data folders other than root.