

**Note:** When graphs are redrawn in live mode, *autoscaling is not done*.

To specify a trace in a graph as being live you must use the **live** keyword with the `ModifyGraph` command. There is no dialog support for this setting.

```
ModifyGraph live(traceName)= mode
```

*Mode* can be 0 or 1. Zero turns live mode off for the given trace.

WaveMetrics provides a demo experiment that generates and displays synthetic data. You should use this experiment to get a feel for the performance you might expect on your particular computer as a function of the window size, number of points in the live wave, and the live modes. To run the demo, choose `File→Example Experiments→Feature Demos→Live Mode`.

Although live mode 1 is not restricted to unity thickness solid lines or dots modes, you will get the best performance if you do use these settings.

### Quick Append

Another feature that may be of use is the quick append mode. It is intended for applications in which a data acquisition task creates new waves periodically. It permits you to add the new waves to a graph very quickly. To invoke a quick append, use the `/Q` flag in an `AppendToGraph` command. There is no dialog support for this setting.

A side effect of quick append is that it marks the wave as not being modified since the last update of graphs and therefore prevents other graphs containing the same wave, if any, from being updated. For a demo, choose `File→Example Experiments→Feature Demos→Quick Append`.

## Graph Preferences

Graph preferences allow you to control what happens when you create a new graph or add new traces to an existing graph. To set preferences, create a graph and set it up to your taste. We call this your *prototype* graph. Then choose `Capture Graph Prefs` from the `Graph` menu.

Preferences are normally in effect only for *manual* operations, not for automatic operations from Igor procedures. This is discussed in more detail in Chapter III-18, **Preferences**.

When you initially install Igor, all preferences are set to the factory defaults. The dialog indicates which preferences you have not changed by displaying “default” next to them.

The `Window Position and Size` preference affects the creation of new graphs only. New graphs will have the same size and position as the prototype graph.

The `Page Setup` preference is somewhat unusual because all graphs share the same page setup settings, as shown in the `Page Setup` dialog. The captured page setup is already in use by all other graphs. The utility of this category is that new *experiments* will use the captured page setup for graphs.

The “`XY Plots:Wave Styles`” preference category refers to the various wave-specific settings in the graph, such as the line type, markers and line size, set with the `Modify Trace Appearance` dialog. This category also includes settings for waveform plots. Each captured wave style is associated with the index of the wave it was captured from. The index of the first wave displayed or appended to a graph is 0, the second appended wave has an index of 1, and so on. These indices are the same as are used in style macros. See **Graph Style Macros** on page II-350.

If preferences are on when a new graph with waves is created or when a wave is appended to an existing graph, the wave style assigned to each is based on its index. The wave with an index of 2 is given the captured style associated with index 2 (the third wave appended to the captured graph).

You might wonder what style is applied to the fifth and sixth waves if only four waves appeared in the graph from which wave style preferences were captured. You have two choices; either the factory default style is used, or the styles repeat with the first wave style and then the second style. You make this choice