

Table Shortcuts

To view table keyboard navigation shortcuts, see **Table Keyboard Navigation** on page II-237.

Action	Shortcut (<i>Macintosh</i>)	Shortcut (<i>Windows</i>)
To add columns for existing waves	Choose Append to Table from the Table menu.	Choose Append to Table from the Table menu.
To create new numeric waves	Click in the first unused column and enter a number or paste numeric data from the clipboard. Copy all or part of an existing numeric wave and paste in the first unused column.	Click in the first unused column and enter a number or paste numeric data from the clipboard. Copy all or part of an existing numeric wave and paste in the first unused column.
To create new text waves	Click in the first unused column and enter nonnumeric text or paste text data from the clipboard. Copy all or part of an existing text wave and paste in the first unused column.	Click in the first unused column and enter nonnumeric text or paste text data from the clipboard. Copy all or part of an existing text wave and paste in the first unused column.
To kill one or more waves	Select the waves in the table and choose Kill Waves from the Table pop-up menu.	Select the waves in the table and choose Kill Waves from the Table pop-up menu.
To rename a wave	Select one or more cells from the wave and choose Rename from the Table pop-up menu.	Select one or more cells from the wave and choose Rename from the Table pop-up menu.
To redimension a wave	Select one or more cells from the wave and choose Redimension from the Table pop-up menu.	Select one or more cells from the wave and choose Redimension from the Table pop-up menu.
To add cells to a column	Click in a row and choose Insert Points from the Table pop-up. Press Command-Shift-V. This does an insert-paste. Click in the insertion cell at the end of the column and enter a value or paste.	Click in a row and choose Insert Points from the Table pop-up. Press Ctrl+Shift+V. This does an insert-paste. Click in the insertion cell at the end of the column and enter a value or paste.
To remove cells from a column	Select the cells to be removed and choose Delete Points from the Table pop-up. Press Command-X (Cut).	Select the cells to be removed and choose Delete Points from the Table pop-up. Press Ctrl+X (Cut).
To identify a particular column	Press Command-Option-Control and click in the column.	Press Shift-F1 to summon context-sensitive help and then click in the column.
To modify column styles	Double-click a column title (goes to Modify Columns dialog). Select one or more columns and click the Table pop-up menu. Press Control and click the column title.	Double-click a column title (goes to Modify Columns dialog). Select one or more columns and click the Table pop-up menu. Right-click the column title.
To modify the style of all columns at once	Press Shift while using the Table pop-up menu. Press Shift-Option to modify all but the Point column.	Press Shift while using the Table pop-up menu. Press Alt+Shift to modify all but the Point column.

Chapter II-12 — Tables

Action	Shortcut (<i>Macintosh</i>)	Shortcut (<i>Windows</i>)
To modify the style of a subset of the data columns of multidimensional wave	Select the columns to modify, press Command, then click the Table menu or the Table pop-up menu, then make a selection.	Select the columns to modify, press Ctrl, then click the Table menu or the Table pop-up menu, then make a selection.
To select entire rows	Drag in the point number column.	Drag in the point number column.
To select entire columns	Drag on the column titles.	Drag on the column titles.
To select all of the columns of a given wave	Press Command and click the column title of a data column from the wave. This selects all data columns from the wave. Press Command and click the column title of an index or label column from a wave to select that column plus all data columns.	Press Ctrl and click the column title of a data column from the wave. This selects all data columns from the wave. Press Ctrl and click the column title of an index or label column from a wave to select that column plus all data columns.
To change a wave's position	Press Option and drag the column title to the new position.	Press Alt and drag the column title to the new position.
To adjust the width of a column	Drag the border at the right of the column title. Press Shift to resize all columns except the point column. Press Command while clicking to set the width of a single data column of a multidimensional wave.	Drag the border at the right of the column title. Press Shift to resize all columns except the point column. Press Ctrl while clicking to set the width of a single data column of a multidimensional wave.
To bring the target cell into view	Click in the cell ID in the top left corner of the table.	Click in the cell ID in the top left corner of the table.